



# MONARK

PRODUCTS FOR LIFE AND PERFORMANCE

NEW  
PRODUCT  
2016!

## Monark LC7TT novo

### How do you compete?

For years athletes and teams have competed in Time Trials against the clock for the fastest time. With the Monark LC7 TT this now becomes a reality in a laboratory environment. Pushing testing and setting up training protocols to form new standards. Whether you want to set the fastest time over a given distance, cover the longest distance in a set time or set the highest workload average, the LC7 TT accurately measure what really matters.



[monarkexercise.se](http://monarkexercise.se)

MONARK LC7 TT novo

# Whats new?



## Workload control

Workload control by +/- buttons placed at the handlebars where the level of each step being selected in the software to create the most individualized settings possible.



## Handle-bars

Easy to change handlebars for road (standard), mtb (option) and triathlon(option), with workload control levers that's fits each one.



## Software module

A specific module for the Monark Test Software for managing the settings and analyzing tools for Time Trial /Open-end testing. The software is developed in cooperation with the Swedish School of Sports and Health Sciences.

## HIGH PERFORMANCE TEST



# MONARK LC7 TT novo

## MAIN FEATURES

- Designed for TimeTrial/ Open-end testing
- Workload control from bars or PC/ECG
- Monark novo touch display
- Monark novo control unit
- Super fast motor
- Perfect fit frame for true cycling position

## THE BASICS

- Perfect cycle feeling
- Easy to use your own saddle, handlebar and pedals
- Speed independent and dependent power control
- Handlebars are easy to change to suit different tests
- Easy to control and calibrate
- Can be controlled by two external devices (RS232 or USB)
- BLE Smart and ANT+ HR systems
- Easy to connect to external devices
- Workload range 4-1400W (200rpm) (8-2800W option)
- Protocols: Full Monark Test Software functionality
- Training mode novo display: Watt, KPM/min, VO<sub>2</sub>, kP and N



Max user weight 180 kg

Medical Device class IIa

Monark Test Software





MONARK EXERCISE AB

PHONE +46 281 594 940 FAX +46 281 719 81 E-MAIL [info@monarkexercise.se](mailto:info@monarkexercise.se) WEBB [monarkexercise.se](http://monarkexercise.se)